



AED
Monthly Inspection form

At a minimum, monthly inspections must verify the following:

Inspection Items	Pass	Fail
Is the unit clean, undamaged, free of excessive wear?		
Are there any cracks or loose parts in the housing?		
Verify electrodes are connected to the AED and sealed in their package. Replace if expired.		
Are all cables free of cracks, cuts, and exposed or broken wires?		
Turn the AED on and off and verify the green check indicated ready for use.		
Batteries within expiration date. Replace if expired.		
Check for adequate supplies.		